DENTASEPT SH PRO

MEDICAL DEVICES

DENTASEPT®

SH PRO

INDICATIONS

Foam detergent disinfectant for the treatment of medical devices in the dental environment.

ADVANTAGES

- ✓ New improved formulation
- ✓ Ready to use solution
- ✓ Formulated without alcohol, compatible with all synthetic materials, plastics etc
- ✓ Without perfume, colourant or allergy promoting substances
- ✓ Foamer spray limits aerosolization of the product
- ✓ Effective in just 1 minute against all enveloped viruses¹
- ✓ Medical Device Class IIb (Directive 93/42/EEC modified)

¹ The Robert Koch Institute recommends the use of disinfectants claiming limited virucidal activity to interrupt the infection-chain. According to the RKI, limited virucidal activity (effective against enveloped viruses) corresponds to an activity against Vaccinia virus and BVDV.

MICROBIOLOGICAL PROPERTIES AT 20°C

Bactericide: EN 1040, EN 13727, EN 13697 - 2 minutes.
MRSA: EN 13727, EN 13697 - 2 minutes.
Yeasticide: EN 1275, EN 13624, EN 13697 - 5 minutes.
Active according to EN 14476 against Rotavirus, Herpes virus, PRV (surrogate of HBV), BVDV (surrogate of HCV), Vaccinia virus - 1 minute.
Active according to EN 14476 against polyomavirus SV40 - 5 minutes.

COMPOSITION

Didecyldimethylammonium chloride, excipients.

PACKAGING

750ml Foamer Spray.
5 litre Container.

DENTASEPT® SH PRO

PROTOCOL FOR USE

Foam detergent disinfectant for the treatment of medical devices in the dental environment.

BEFORE USE

Wear appropriate gloves.


HOW TO USE

- 1 Spray Dentasept® SH PRO on to a non-woven cloth or directly on to the area to be cleaned and disinfected.
- 2 Evenly distribute the product (+/- 6 ml/m² or approx. 10 sprays) with a non-woven cloth to ensure that the complete area is treated.
- 3 Observe the contact time specified to achieve the level of activity required.
Rinsing is not necessary unless the medical device comes into contact with skin or mucous membranes.

PRECAUTIONS

Dangerous - respect the precautions for use. Medical Device Class IIb, CE 0459. Before use, read the label and information concerning the product.

Store between +5°C and +35°C.

Summary

Regulation and General Information Brochure	Page	7
Composition	Page	9
Stability and storage conditions	Page	11
Studies and expertises		
- Microbiology	Page	13
- Corrosion	Page	17
- Compatibilities	Page	19
Regulation	Page	21

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Regulation and general information

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***DENTASEPT SH PRO** is designed, produced and controlled by the Laboratoires ANIOS, certified by the AFAQ Organism under the number 1995/3723 in accordance with the ISO 9001 Quality System*

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***DENTASEPT SH PRO** has the CE mark for medical devices in the category class IIb, in accordance with the Directive 93/42/CE*

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***DENTASEPT SH PRO** formulation fulfils the requirements of Regulation (EC) No 648/2004 on detergents and its amendments*

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***DENTASEPT SH PRO** is labelled in accordance with the European regulation related to the classification and labelling of chemical products*

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***DENTASEPT SH PRO** meets our voluntary eco-design commitment included in our ANIOSAFE charter*

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Pure product ready to use :

- *Storage between +5°C and +35°C.*
- *Stability:*
- *Non opened bottle: 2 years from the manufacturing date indicated on the label.*
- *Spray: preserving of the 2 years if the product is stored in a container closed by the sprayer.*

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It is shown that, according to the NF S 94-402-1 standard (soaking method), DENTASEPT SH PRO does not corrode stainless steel 304 L in the conditions of use recommended by the Laboratoires ANIOS.

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*Red label (on the left): FORMULA 2419
White label (on the right): demineralized water*

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Not concerned by biodegradability notion

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0,5 %¹

*Readily biodegradable organic material*²

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Biodegradability within 28 days or less:

- ≥ 60% according to respirometric methods³ (assessment of CO₂ production or oxygen consumption)
- ≥ 70% according to methods measuring dissolved organic carbon³ (DOC) (assessment of DOC die away)

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Substances which not fulfil the readily biodegradability criteria (but possibly displaying inherent biodegradability criteria⁴), or substances without any biodegradability data.

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DENTASEPT SH PRO contains 100.0 % of inorganic material and of readily biodegradable organic material.

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HAZARDS IDENTIFICATION

Health

*This mixture does not present a health hazard with the exception of possible occupational exposure thresholds**.**

Environment

*Hazardous to the aquatic environment - Chronic hazard, Category 3 (Aquatic Chronic 3, H412).**

** Classification drawn up in accordance with European regulation related to the classification and labelling of chemical products.*

*** please refer to SDS §8 for TWA (Time-Weighted Average)/STEL (Short-term exposure limit) of concerned substances present in the product.*

PERSONAL PROTECTION UNDER THE CONDITIONS OF THE INTENDED USE

Refers exclusively to this product. If the product is handled and/or exposed simultaneously with other chemical agents, these must be taken into consideration when choosing personal safety equipment.

Respiratory protection

Not applicable in normal conditions of use.

Hand protection

When handling this product, wear suitable gloves.

Nitrile, latex or vinyl gloves.

Gloves must be replaced immediately if they show signs of wear and tear.

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The Safety Data Sheet (SDS) of this product is available on our Web site www.anios.com. In addition to visualization and/or impression, it is possible for you to download it.

The consultation of this SDS via internet enables you to receive later on by email any updated version of this document.

In the absence of an Internet access, the SDS remains available on request near our laboratories (tel.: +33 3 20 67 67 67 - fax: +33 3 20 67 67 68)

* an access is also possible at the site <http://www.quickfds.com>